

Bacillo[®] 30 Foam

Alcohol-based
rapid disinfection



Aldehyde-free, material-friendly rapid disinfectant for disinfecting sensitive surfaces and medical equipment.

Bacillo[®] 30 Foam

Product properties

- ready-to-use disinfectant solution
- acts rapidly
- especially material-friendly, also suitable for sensitive materials such as Makrolon[®], Plexiglas[®] and Polysulfon
- excellent wetting
- handy spray bottle with foam spray head ensures user-friendly application
- colourant and fragrance-free

Composition

Active ingredients: Ethanol 140 mg/g; propan-2-ol 100 mg/g; propan-1-ol 60 mg/g, n-alkyl-aminopropyl-glycine (CAS number 1397 34-65-9) 5 mg/g

Microbiology

- bactericidal
- yeasticidal
- virucidal against enveloped viruses (incl. HBV, HIV, HCV)
- MNV
- polyoma- and rotaviruses

Areas of application

Bacillo[®] 30 Foam is suitable for the rapid spray/wipe disinfection of sensitive washable surfaces, e.g.:

- medical equipment and devices that come under the Medical Device Directive (MDD)

Surfaces acc. to Biocidal Products Directive (BPD), e.g.:

- sensitive synthetic surfaces such as synthetic leather, Makrolon[®], acrylic glass and polysulfone
- working surfaces in clinics, primary health care, dental surgeries and ambulances
- areas that require rapid effect, e.g. production lines in order to reduce down times
- areas difficult to access, such as threads, stop valves, etc.
- displays/ monitors and keyboards/ control panels
- solaria

Directions for use

Bacillo[®] 30 Foam is ready to use. Sufficiently wet the surface to be disinfected and wipe it with a clean disposable cloth. Keep the surface moist during the entire exposure time.

For an even application use the especially developed foam spray head, which, in addition, facilitates the optical control of the coverage.

For the disinfection of vertical surfaces, e.g. displays and monitors, and uneven surfaces such as keyboards and control panels, we recommend spraying a sufficient amount of Bacillo[®] 30 Foam on a clean disposable cloth to completely and thoroughly wipe the surfaces to be disinfected.

Please make sure that no disinfection solution gets inside the device. Where applicable, read the operating manual of the device to be disinfected for detailed information. Do not use for the disinfection of invasive medical devices.

The amount of disinfection solution applied must not exceed 50 ml per m² of surface to be treated. The total amount applied per room must not exceed 100 ml per m² of room floor space.

Do not allow product to reach water systems undiluted.

Use disinfectants safely. Always read the label and product information before use.

Dosage

Bacteria/Fungi

Efficacy according to EN Phase 2/Step 1 (suspension tests), tested under clean / dirty conditions	
Bactericidal (EN 13727)	
- dirty conditions	30 sec.
Yeasticidal (EN 13624)	
- dirty conditions	30 sec.
Certified application recommendations for prophylactic wipe disinfection from the Association for Applied Hygiene (VAH). Based on suspension and practical tests, tested under clean conditions (i.e. optically clean surfaces) / dirty conditions (i.e. visibly contaminated surfaces)	
Bactericidal/Yeasticidal	
- dirty conditions	5 min.
Rapid disinfection (in accordance with the German Society for Hygiene and Microbiology [DGHM]); Based on suspension and practical tests; tested under clean / dirty conditions	
Bactericidal/yeasticidal	
- dirty conditions	30 sec.

Viruses

Efficacy against viruses (German Society for the Control of Viral Diseases [DVV])	
Virucidal against enveloped viruses (incl. HBV, HIV, HCV)	30 sec.
Appraised efficacy against non-enveloped viruses (DVV)	
Polyomavirus	5 min.
Appraised efficacy against non-enveloped viruses (in accordance with EN)	
MNV	
- clean conditions	5 min.
- dirty conditons	10 min.
Rotavirus	1 min.

Listing

- Certification/list issued by the Association for Applied Hygiene (VAH) (application submitted)
- CE-labelling in accordance with the Medical Device Directive (MDD)



Bacillo[®] 30 Foam

Chemical-physical data

Flashpoint (DIN 51755):	30 °C
Density (20 °C):	approx. 0.96 g/cm ³
pH-value (20 °C):	approx. 8
Refraction index n_D^{20} :	approx. 1.36

Presentation

750 ml bottle*, 5 litre canister and 200 litre barrel

* comes with 2 foam spray heads per sales unit of 8 bottles

Additional application aids on demand.

Labelling in acc. with the 1999/45/EG guideline: Xi Irritant. Flammable. Irritating to eyes. Vapours may cause drowsiness and dizziness. For correct use, the basic measures of occupational health and safety corresponding to protection level 2, Art. 9 of the Ordinance on Hazardous Materials have to be observed. At temperatures above its flashpoint of 30 °C, Bacillo[®] 30 Foam may form explosive atmospheres. Do not breathe vapour/aerosol. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Use only in well-ventilated areas.

02-11

