

Safety Data Sheet

according to 1907/2006/EC, Article 31

Date of issue: 10.12.2008

Revision: 10.12.2008

1 Identification of substance:

- **Product details:**
- **Trade name:** Bodedex® forte
- **SAP-Code:** EN R10214
- **Application of the substance / the preparation** Instrument cleaner
- **Manufacturer/Supplier:**
BODE Chemie GmbH & Co. KG
Melanchthonstraße 27
D-22525 Hamburg (Germany)
Tel.: +49 (0)40 / 54 00 60
- **Informing department:** .
- **Emergency information:**
Giftnotruf Göttingen (Germany)
24h-Tel.: +49 (0)551 / 1 92 40

2 Hazards identification

- **Hazard designation:**



Xi Irritant

- **Information pertaining to particular dangers for man and environment**
R 38 Irritating to skin.
R 41 Risk of serious damage to eyes.
- **Classification system**
The classification is in line with current EC lists. It is expanded, however, by information from technical literature and by information furnished by supplier companies.

3 Composition/information on ingredients

- **Chemical characterization**
- **Description:** Mixture of the substances listed below with harmless additions.

- **Dangerous components:**

nonionic surfactant	☒ Xn, ☒ Xi; R 22-41	5-10%
amphoteric surfactant	☒ Xi, ☒ N; R 36/38-50	1-5%

- **Regulation (EC) No 648/2004 on detergents / Labelling for contents**

non-ionic surfactants	5 - 15%
phosphonates, amphoteric surfactants	< 5%
(ethylenedioxy)dimethanol, Mixture of 5-Chloro-2-methylisothiazol-3(2H)-one and 2-methylisothiazol-3(2H)-one with magnesium chloride and magnesium nitrate	

- **Additional information** For the wording of the listed risk phrases refer to section 16.

4 First aid measures

- **After skin contact**
Instantly rinse with water.
If skin irritation continues, consult a doctor.
- **After eye contact**
Rinse opened eye for several minutes under running water.
Call a doctor immediately.

(Contd. on page 2)

Safety Data Sheet

according to 1907/2006/EC, Article 31

Date of issue: 10.12.2008

Revision: 10.12.2008

Trade name: Bodedex® forte

(Contd. of page 1)

- **After swallowing**
Rinse out mouth and then drink plenty of water.
In case of persistent symptoms consult doctor.

5 Fire fighting measures

- **Suitable extinguishing agents** Use fire fighting measures that suit the environment.
- **Protective equipment:** No special measures required.

6 Accidental release measures

- **Person-related safety precautions:** Wear protective clothing.
- **Measures for environmental protection:** Do not allow to enter drainage system, surface or ground water.
- **Measures for cleaning/collecting:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

7 Handling and storage

- **Handling**
- **Information for safe handling:** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Storage**
- **Requirements to be met by storerooms and containers:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep container tightly sealed.

* 8 Exposure controls and personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Components with critical values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the compilation were used as basis.
- **Personal protective equipment**
- **General protective and hygienic measures**
Keep away from foodstuffs, beverages and food.
Instantly remove any soiled and impregnated garments.
Wash hands during breaks and at the end of the work.
- **Breathing equipment:** Not required.
- **Protection of hands:**
Impervious gloves.
Only use chemical-protective gloves with CE-labelling of category 3.
Check protective gloves prior to each use for their proper condition.
After use of gloves apply skin-cleaning agents and skin cosmetics.
- **Material of gloves**
 1. Nitrile rubber, NBR
 2. Butyl rubber, BR
Recommended thickness of the material: 1. 0,4 mm / 2. 0,7 mm
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 3)

Safety Data Sheet

according to 1907/2006/EC, Article 31

Date of issue: 10.12.2008

Revision: 10.12.2008

Trade name: **Bodedex® forte**

(Contd. of page 2)

- **As protection from splashes gloves made of the following materials are suitable:**
 1. Nitrile rubber, NBR
 2. Butyl rubber, BR
- **Not suitable are gloves made of the following materials:**
 - Natural rubber, NR
 - PVA gloves
- **Eye protection:** Tightly sealed safety glasses.

9 Physical and chemical properties:

· General Information

Form:	fluid
Colour:	Light yellow
Smell:	mild

· Change in condition

Melting point/Melting range: not determined
Boiling point/Boiling range: 100°C

· **Flash point:** > 100°C

· **Self-inflammability:** product is not selfigniting

· **Danger of explosion:** product is not explosive

· **Density at 20°C** 1,05 g/cm³

· **Solubility in / Miscibility with Water:** fully miscible

· **pH-value at 20°C:** 8

10 Stability and reactivity

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Dangerous reactions** No dangerous reactions known.
- **Dangerous products of decomposition:** No dangerous decomposition products known.

11 Toxicological information

- **Acute toxicity:**
- **Primary irritant effect:**
 - **on the skin:** Irritant for skin and mucous membranes.
 - **on the eye:** Strong irritant with the danger of severe eye injury.
- **Sensitization:** No sensitizing effect known.
- **Additional toxicological information:**
 The product shows the following dangers according to the calculation method of the General EC Classification Guidelines for Preparations as issued in the latest version:
 Irritant

12 Ecological information:

- **Information about elimination (persistence and degradability):**
- **Other information:**
 The surfactant(s) contained in this preparation complies(comply) with the biodegradability criteria as laid down in Regulation (EC) No.648/2004 on detergents. Data to support this assertion are held at the disposal of the competent authorities of the Member States and will be made available to them, at their direct request

(Contd. on page 4)

Safety Data Sheet

according to 1907/2006/EC, Article 31

Date of issue: 10.12.2008

Revision: 10.12.2008

Trade name: Bodedex® forte

(Contd. of page 3)

or at the request of a detergent manufacturer.

· **Additional ecological information:** -

13 Disposal considerations

· **Product:**

· **Recommendation**

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.
Must be specially treated under adherence to official regulations.

· **Uncleaned packagings:**

· **Recommendation:** Dispose of packaging according to regulations on the disposal of packagings.

· **Recommended cleaning agent:** Water, if necessary with cleaning agent.

14 Transport information

· **Land transport ADR/RID and GGVS/GGVE (cross-border/domestic)**

· **ADR/RID-GGVS/E Class:** -

· **Maritime transport IMDG/GGVSea:**

· **IMDG/GGVSea Class:** -

· **Air transport ICAO-TI and IATA-DGR:**

· **ICAO/IATA Class:** -

· **Transport/Additional information:** Not dangerous according to the above specifications.

15 Regulatory information

· **Designation according to EC guidelines:**

The product has been classified and labelled in accordance with EC Directives / Ordinance on Hazardous Materials (GefStoffV).

· **Code letter and hazard designation of product:**



Xi Irritant

· **Risk phrases:**

38 Irritating to skin.

41 Risk of serious damage to eyes.

· **Safety phrases:**

26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

37/39 Wear suitable gloves and eye/face protection.

16 Other information:

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant R-phrases**

22 Harmful if swallowed.

36/38 Irritating to eyes and skin.

41 Risk of serious damage to eyes.

50 Very toxic to aquatic organisms.

· *** Data compared to the previous version altered.**